

## Single Use Manual Resuscitator



Effective. Reliable. Robust.

**flexicare**  
TOTAL QUALITY - TOTAL CARE

Flexicare's range of ergonomically designed Single Use Manual Resuscitators with higher levels of Oxygen delivery and pre-fitted components in quick-to-open packaging, provide fast effective care, even in the most difficult situations.



Ventilation

Option of 60 or 40 cmH<sub>2</sub>O Pressure Relief Valves to prevent lung over-inflation and minimise gastric insufflation.

Highest oxygen concentrations at volumes recommended in the ERC 2015 Resuscitation Guidelines<sup>1</sup> and International Standards<sup>2</sup>.


| Delivered Oxygen Concentration % |         |         |         |         |          |
|----------------------------------|---------|---------|---------|---------|----------|
| Oxygen Flow Rates                | 2 L/min | 4 L/min | 6 L/min | 8 L/min | 15 L/min |
| Adult                            | 40%     | 56.5%   | 80%     | 98%     | 99.5%    |
| Paediatric                       | 45%     | 70%     | 86%     | 98%     | 100%     |
| Infant                           | 67%     | 85%     | 91%     | 98%     | 100%     |

PACKAGING


Quick-to-open packaging and pre-fitted mask minimise set up time. Additional devices, such as an adjustable PEEP valve or Manometer, are fitted without removing the Resuscitator from the patient, so ventilation can begin immediately and continue uninterrupted.



The carry bag handle can be used to hang the Resuscitator in an easily accessible location, taking up minimal valuable storage space.



Clear packaging is large enough to expand the device for routine testing, whilst the tear-off strip assures the user that the contents are uncontaminated.



Once used, the contaminated Resuscitator is resealed in the carry bag for safe transport and disposal.

Capnography

ERC 2015 Resuscitation Guidelines<sup>1</sup> states that "Waveform capnography must be used for confirming tracheal tube placement and monitoring ventilation rate" (In-Hospital CPR) and "must be used whenever tracheal intubation is undertaken" (Prehospital CPR). The ERC recommend capnography during all CPR to monitor the quality of chest compressions, avoid hyperventilation, identify ROSC and prognostication, and "considered for patients ventilated by an SGA".





## USABILITY

High oxygen concentrations at optimal tidal volumes, with the option of a blue 60 cmH<sub>2</sub>O or yellow 40 cmH<sub>2</sub>O Pressure Relief Valve, limits excessive ventilation, safeguarding the patient from harm.

### Lightweight Construction

Low resistance valves and fast recoil bag reduce user fatigue, even during lengthy ventilations.

### Grip & Control

Textured surface and optional strap gives grip and control of the ambidextrous bag, even when wearing gloves and in wet conditions.

### Optional Strap

Facilitating single handed bag use, the strap frees the other hand to secure the mask or airway device. Double handed compression for adults remains possible if required.

### Colour Coded Pressure Relief Valve

Located centrally, the Pressure Relief Valve can be operated by either hand and the setting is visible from both sides and above.

### Versatility

All fittings are International Standard<sup>3,4</sup> mask, ET Tube, Laryngeal Mask and Oxygen line connections, avoiding potential delays.


### Double Swivel Patient End Connector

The Valve Assembly rotates 360° in two axes for an effective mask seal and bag position, whatever the approach angle.

## SAFETY

### Clear Valve Assembly

Unrestricted view of correct valve function, and reduced bulk so the user can verify mask or airway position and observe the patient.



### Non-Tangle Oxygen Tubing

Star lumen Oxygen tubing maintains gas flow even if kinked or crushed. Permanently bonded to avoid accidental disconnection.



### Textured Reservoir Bag

The non-stick textured inner surface ensures complete inflation and correct bag function.



### Optional CO<sub>2</sub> Monitoring

Located close to the patient airway for maximum capnography trace accuracy & reliability, the capless self-sealing port eliminates gas leak risk when not in use.

### Optional MDI Port

The MDI port is angled so medication is directed into the patient airway, maximising effectiveness and reducing waste.

### Optional Manometer

Available in pack or as a separate item, the single use manometer gives a breath-by-breath read-out of delivered pressures.

## ORDERING INFORMATION



Adult



Paediatric



Infant

| PART NUMBER | DESCRIPTION                                 | + PRESSURE RELIEF VALVE | + OPTIONAL PORTS      | + ACCESSORY | MASK SIZE | QUANTITY |
|-------------|---|-------------------------|-----------------------|-------------|-----------|----------|
| 038-71-940  | Adult Single Use Resuscitator               | 60cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | -           | 5         | 5        |
| 038-71-941  | Adult Single Use Resuscitator               | 60cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | Manometer   | 5         | 5        |
| 038-71-942  | Adult Single Use Resuscitator               | 60cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | PEEP Valve  | 5         | 5        |
| 038-71-943  | Adult Single Use Resuscitator with Strap    | -                       | -                     | -           | 5         | 5        |
| 038-71-944  | Adult Single Use Resuscitator without Strap | -                       | -                     | -           | 5         | 5        |
| 038-71-945  | Adult Single Use Resuscitator               | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | -           | 5         | 5        |
| 038-72-940  | Paediatric Single Use Resuscitator          | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | -           | 2 + 1     | 5        |
| 038-72-941  | Paediatric Single Use Resuscitator          | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | Manometer   | 2 + 1     | 5        |
| 038-72-942  | Paediatric Single Use Resuscitator          | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | PEEP Valve  | 2 + 1     | 5        |
| 038-72-946  | Paediatric Single Use Resuscitator          | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | -           | 3         | 5        |
| 038-73-940  | Infant Single Use Resuscitator              | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | -           | 1 + 0     | 5        |
| 038-73-941  | Infant Single Use Resuscitator              | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | Manometer   | 1 + 0     | 5        |
| 038-73-942  | Infant Single Use Resuscitator              | 40cmH <sub>2</sub> O    | CO <sub>2</sub> & MDI | PEEP Valve  | 1 + 0     | 5        |
| 038-71-966  | Single Use Manometer                        |                         |                       |             |           | 5        |

### – Notes –

<sup>1</sup> European Resuscitation Council Guidelines for Resuscitation 2015. Resuscitation, Vol. 95. October 2015.

<sup>2</sup> BS EN ISO 10651-4:2009 Lung Ventilators. Part 4: Particular requirements for operator-powered resuscitators.

<sup>3</sup> ISO 5356-1 Anaesthetic and respiratory equipment - Conical connectors - Part 1: Cones and Sockets

<sup>4</sup> ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment - Part 1: General Requirements



www.flexicare.com

FC3400 | 022017

#### Flexicare UK

Flexicare Medical Limited  
Cynon Valley Business Park,  
Mountain Ash, CF45 4ER, UK  
☎ +44 (0)1443 474647  
✉ +44 (0)1443 474222  
@ enquiries@flexicare.com

#### Flexicare India

Flexicare Medical India Pvt.Ltd  
North Phase, Developed Plot  
No.10A, T.V.K. Industrial Estate,  
Ekkattuthangal,  
Chennai-600 097, India  
☎ +91 44 4554 0220  
@ enquiries@flexicare.com

#### Flexicare USA

Flexicare Inc  
15281 Barranca Pkwy, Unit D,  
Irvine, CA 92618, USA  
☎ 949-450 9999  
✉ 949-450 9992  
@ inquiries@flexicare.com

#### Flexicare China

Flexicare Medical Dongguan Limited  
No. B-15, Xicheng Industrial  
Zone I, Hengli Town, Dongguan City,  
Guangdong Province  
China 523460  
☎ +86 769 83713371  
@ enquiries@flexicareasia.com

#### Flexicare MEA

Flexicare Middle East & Africa  
Um Uthaina, Shat Al Arab Street,  
Bld. 9, 4th Floor, Amman, Jordan  
☎ +962 6554 8944  
✉ +962 6554 8945  
@ enquiries@flexicare.com

#### Flexicare Australia

Flexicare Australia Pty Ltd  
Suite 29A, Level 2,  
541 Blackburn Rd  
Mt Waverley, Victoria 3149  
Australia  
☎ +61 450 977 007  
@ enquiries@flexicare.com